



Acr-Bis (29:1) Datasheet

Introduction:

30% Acr-Bis (29:1) is the aqueous solution contains 30% acrylamide-bisacrylamide, and the ratio of acrylamide and bisacrylamide is 29:1.

30% Acr-Bis (29:1) Commonly used in the preparation of acrylamide gel (PAGE gel), including SDS-PAGE gel, etc., can be used for protein or nucleic acid separation. It is a common reserve liquid in laboratory.

Contents:

Cat No.: BD0036

Acr-Bis(29:1) 100ml

Attention:

30% Acr-Bis(29:1) is toxic. Please pay attention to proper protection when using it.

The product is limited to the scientific research of professionals. It can not be used for clinical diagnosis or treatment. It can not be used for food or medicine. It can not be stored in ordinary houses.

For your safety and health, please wear test clothes and disposable gloves.

Storage & Shelf life:

Store at 4°C. ,Period of validity for one years .

Bioworld Technology, Inc.

1660 South Highway 100, Suite 500 St. Louis Park, MN
55416,USA. Email: info@bioworld.com
Tel: 6123263284 Fax: 6122933841

Bioworld technology, co, Ltd.

No 9, weidi road Qixia District Nanjing, 210046,
P, R.China. Email: info@biogot.com
Tel: 0086-025-68037686 Fax:0086-025-68037575