



Tris- Hydrochloric Acid Spacer Gel Buffer Datasheet

Introduction:

0.5mol/L Tris-Hcl, PH6.8 used to prepare SDS-PAGE or native PAGE superstratum gel.

Contents:

Cat No.: BD0037-2

Tris- Hydrochloric Acid Spacer Gel Buffer 60ml

Attention:

The product is limited to the scientific research of professionals. It can not be used for clinical diagnosis or treatment. It can not be used for food or medicine. It can not be stored in ordinary houses.

For your safety and health, please wear test clothes and disposable gloves.

Storage & Shelf life:

Store at 4°C. ,Period of validity for six months .

Bioworld Technology, Inc.

1660 South Highway 100, Suite 500 St. Louis Park, MN

55416,USA.

Email: info@bioworld.com

Tel: 6123263284

Fax: 6122933841

Bioworld technology, co, Ltd.

No 9, weidi road Qixia District Nanjing, 210046,

P, R.China.

Email: info@biogot.com

Tel: 0086-025-68037686 Fax:0086-025-68037575