



Tris- Glycine-SDS Electrophoretic Buffer 10x

Introduction:

Tris- Glycine Electrophoretic Buffer (SDS-PAGE Electrophoretic Buffer). This reagent is a 10-fold concentrated solution, which needs to be diluted by distilled water or deionized water for 10-fold use.

Contents:

Cat No.: BD0061

Tris- glycine -SDS Electrophoretic Buffer10x 1L

Procedure:

Diluted distilled water or deionized water diluted to 1 X before use.

Attention:

1. Occasionally, there are turbidity or precipitation, which can be used after heating and dissolving.
2. If the amount of each use is very small, it can be properly sub-packed before use.
3. For your safety and health, please wear test clothes and disposable gloves.

Storage & Shelf life:

Store at Room temperature ,Period of validity for three months.

Bioworld Technology, Inc.

1660 South Highway 100, Suite 500 St. Louis Park, MN

55416,USA.

Email: info@bioworld.com

Tel: 6123263284

Fax: 6122933841

Bioworld technology, co, Ltd.

No 9, weidi road Qixia District Nanjing, 210046,

P, R.China.

Email: info@biogot.com

Tel: 0086-025-68037686 Fax:0086-025-68037575

For Research Use Only