

Protease inhibitor cocktail (100X)

Introduction:

The product is a general protease inhibitor cocktail that can be used for protein extraction from cells or tissues. It can effectively inhibit various proteases in conventional cell or tissue extracts, such as serine proteases, aminopeptidases, cysteine proteases, threonine and aspartic proteases.

Contents:

Cat No.: BD0075 Protease inhibitor cocktail (100X) : 1ml

Procedure:

The product should be added to lysis buffer at a ratio of 1:100 for use. Apply to lysis buffer before use.

Attention:

The product is limited to the scientific research of professionals. It can not be used for clinical diagnosis or treatment. It can not be used for food or medicine. It can not be stored in ordinary houses.

For your safety and health, please wear test clothes and disposable gloves.

Storage & Shelf life:

Store at -20°C

Bioworld Technology, Inc. 1660 South Highway 100, Suite 500 St. Louis Park, MN 55416,USA. Email: info@bioworlde.com Tel: 6123263284 Fax: 6122933841 Bioworld technology, co, Ltd. No 9, weidi road Qixia District Nanjing, 210046, P, R.China. Email: <u>info@biogot.com</u> Tel: 0086-025-68037686 Fax:0086-025-68037575

For Research Use Only