

## PRODUCT DATA SHEET

---



# Protease inhibitor (PMSF)

### Introduction:

10mM/L PMSF, Working concentration 0.1-1mM

### Contents:

Cat No.: BD0033

Protease inhibitor(PMSF) 2ml

### Procedure:

PMSF is unstable in aqueous solution. Apply to lysis buffer before use.

### Attention:

Toxic, Attention protection.

The product is limited to the scientific research of professionals. It can not be used for clinical diagnosis or treatment. It can not be used for food or medicine. It can not be stored in ordinary houses.

For your safety and health, please wear test clothes and disposable gloves.

### Storage & Shelf life:

Store at -20°C

---

Bioworld Technology, Inc.

1660 South Highway 100, Suite 500 St. Louis Park, MN

55416,USA.

Email: [info@bioworld.com](mailto:info@bioworld.com)

Tel: 6123263284

Fax: 6122933841

Bioworld technology, co, Ltd.

No 9, weidi road Qixia District Nanjing, 210046,

P, R.China.

Email: [info@biogot.com](mailto:info@biogot.com)

Tel: 0086-025-68037686 Fax:0086-025-68037575

For Research Use Only