

PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



5% BSA Blocking Buffer (TBS Solute)

Cat No.: BD5020

Introduction

This product contains Bovine Serum Albumin, TBS, Preservatives and indicator, without Triton X-100. It has been filtered through out 0.22um membrane. This product is suitable for WB, ELISA and immunohistochemistry (or immunofluorescence) blocking. If the buffer color and smell is changed, you'd better to discard it.

Reagents

5% BSA Blocking Buffer (TBS Solute) 100ml

Application

For blocking in IHC, ICC, IF and WB.

Storage & Shelf life

Store at 2-8 °C. Each component is stable for up to 12 Months.

Research Use

For research use only, not for use in diagnostic procedures.

MADE IN CHINA

Bioworld Technology, Inc.
1660 South Highway 100, Suite 500 St. Louis Park, MN
55416, USA. Email: info@bioworld.com
Tel: 6123263284 Fax: 6122933841

Bioworld technology, co, Ltd.
No 9, weidi road Qixia District Nanjing, 210046,
P, R.China. Email: info@biogot.com
Tel: 0086-025-68037686 Fax: 0086-025-68035151