

PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



Frozen Section Protect Buffer

Cat No.: BD5034

Introduction

This product is used for protect frozen section at -20°C , especially for drift section of central nervous system sections. It contains sucrose, PBS, ethylene glycol and preservatives. This protect buffer can prevent the section generate ice crystal that is harmful for antigen.

Reagents

Frozen Section Protect Buffer 250ml

Application

Frozen Section Protect in IHC-Fr.

Storage & Shelf life

Store at $2-8^{\circ}\text{C}$. Each component is stable for up to 6 months.

NOTE

1. The Section Protect Buffer can be frozen at low temperature blew -25°C ;
2. The fresh frozen section has the best antigenicity;
3. Avoid freeze/ thaw cycles;
4. The section should not be store in Protect Buffer over 3 months.

Research Use

For research use only, not for use in diagnostic procedures.

MADE IN CHINA

Bioworld Technology, Inc.
1660 South Highway 100, Suite 500 St. Louis Park, MN
55416, USA. Email: info@bioworld.com
Tel: 6123263284 Fax: 6122933841

Bioworld technology, co, Ltd.
No 9, weidi road Qixia District Nanjing, 210046,
P, R.China. Email: info@biogot.com
Tel: 0086-025-68037686 Fax: 0086-025-68035151