

## PRODUCT DATA SHEET

---



### Wash buffer for WB (25X)

#### Introduction:

The product can be used for washing after incubation of primary or secondary antibodies at Western. Proper washing can reduce the background and enhance the signal-to-noise ratio.

#### Contents:

Cat No.: WBA2400-R6  
25\*TBST Rinse Buffer 1L

#### Procedure:

When the product is in use, it needs to be equipped with 10 \* TBST rinse buffer, and then with 10 \* TBST rinse buffer to be equipped with 1 \* TBST rinse buffer. Taking the make 500 ml 10 \* TBST rinse buffer as an example, First, 200 ml 25 \* rinse buffer, adding 300 ml distilled water, mixing, that is, 10 \* TBST rinse buffer, then 100 ml 10 \* TBST rinse buffer, adding 900 ml distilled water, mixing, that is, 1 \* TBST rinse buffer. 1 \* TBST, can be used for Western washing. The unused solution can be stored at 4°C, and can normally be used within 1 - 2 weeks. If abnormal phenomena such as turbidity and sediment are found in the detergent, it should be discarded.

#### Attention:

The 25 \* TBST rinse buffer is prepared into 1 \* TBST rinse buffer, which should be stored at 4°C and used up within 1 - 2 months. If stored at room temperature, it should be used up within 1-2 weeks. Unconfigured 25\*TBST rinse buffer or 10\*TBST rinse buffer can be stored at room temperature. The product is limited to the scientific research of professionals. It can not be used for clinical diagnosis or treatment. It can not be used for food or medicine. It can not be stored in ordinary houses.

For your safety and health, please wear test clothes and disposable gloves.

#### Storage & Shelf life:

Storage at 4°C is valid for one year. It can be preserved for a long time at -20°C, and can be kept at room temperature for one month.

---

**Bioworld Technology, Inc.**

1660 South Highway 100, Suite 500 St. Louis Park, MN  
55416,USA. Email: [info@bioworld.com](mailto:info@bioworld.com)  
Tel: 6123263284 Fax: 6122933841

**Bioworld technology, co, Ltd.**

No 9, weidi road Qixia District Nanjing, 210046,  
P, R.China. Email: [info@biogot.com](mailto:info@biogot.com)  
Tel: 0086-025-68037686 Fax:0086-025-68037575

**For Research Use Only**